
KORLEIS FÅ RIKTIGE ANALYSERESULTAT PÅ BLODGASSANALYSENE?

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Kjelder

1. Higgins C. Useful tips to avoid preanalytical errors in blood gas testing: pH, pCO₂ and pO₂. www.acutecaretesting.org May 2016
2. Wennecke G. Useful tips to avoid preanalytical errors in blood gas testing: electrolytes. www.acutecaretesting.org Oct 2003
3. Wennecke G. Useful tips to avoid preanalytical errors in blood gas testing: metabolites. www.acutecaretesting.org Jan 2004
4. Plebani M et al. Performance criteria and quality indicators for the pre-analytical phase. *Clin Chem Lab Med* 2015; 53: 943-48
5. CLSI. Procedures for the collection of arterial blood specimens; Approved standard – Fourth Edition. CLSI document GP43-A4. Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19807 USA, 2004.
6. Hajiseyedjavady H *et al.* Less painful arterial blood gas sampling using jet injection of 2% lidocaine: a randomised controlled clinical trial. *Am J Emerg Med* 2012; 30: 1100-04.
7. Weinrich U *et al.* Time to steady state after changes in FiO₂ in patients with COPD. *J Obstructive Pulmon Dis* 2013; 10: 405-10.
8. Higgins C. Central venous blood gas analysis. www.acutecaretesting.org July 2011.
9. Byrne A *et al.* Peripheral venous and arterial blood gas analysis in adults: are they comparable? A systematic review and meta-analysis. *Respirology* 2014; 19: 168-75.
10. Higgins C. Capillary blood gases: to arterialize or not. www.acutecaretesting.org July 2008.
11. Zavorsky G *et al.* Arterial versus capillary blood gases: a meta-analysis. *Respiratory Physiol and Neurobiol* 2007; 155: 268-79.
12. Toffaletti J *et al.* Effect of small air bubbles on changes in pO₂ and blood gas parameters: calculated vs. measured effects. www.acutecaretesting.org July 2012.
13. Dean J *et al.* Effects of syringe material and temperature and duration of storage on the stability of equine arterial blood gas variables. *Veterinary Anaesthesia and Analgesia* 2004; 31: 250-57.
14. Knowles T *et al.* Effects of syringe material sample storage time and temperature on blood gases and oxygen saturation in arterialized human blood samples. *Respiratory Care* 2006; 51: 732-36.
15. Wiwanitikit V. Glass syringes are better than plastic for preserving arterial blood gas for oxygen partial pressure determination: an explanation based on nanomaterial composition. *Int J Nanomedicine* 2006; 1: 223-24.
16. Higgins C. The use of heparin in preparing samples for blood gas analysis. www.acutecaretesting.org April 2007.
17. Hutchinson A *et al.* Too much heparin: possible source of error in blood gas analysis *BMJ* 1983; 287: 1131-32.
18. Chapola V *et al.* Is liquid heparin comparable to dry balanced heparin for blood gas sampling in intensive care unit. *Indian J Crit Care Med* 2014; 18: 14-20.
19. Grenache D *et al.* Integrated and automatic mixing of whole blood: an evaluation of a novel blood gas analyser. *Clin Chim Acta* 2007; 375: 153-57.
20. Dikmen Z *et al.* Specimen rejection in laboratory medicine: Necessary for patient safety? *Biochimica Medica* 2015; 25: 377-85.
21. Astles J *et al.* Pneumatic transport exacerbates interference of room air contamination in blood gas samples. *Arch Pathol Lab Med* 1993; 120: 642-47.
22. Higgins C. Pneumatic tube transport of blood samples – an update. www.acutecaretesting.org Oct 2015.
23. Blonshine S, Alberti R, Olesinski RL. Procedures for the collection of arterial specimens. Approved standard. 3rd ed. NCCLS Publication H11-A3. Wayne, Pa.: NCCLS, 1999; 19, 8.
24. Siggaard-Andersen O, Thode J, Wandrup J. The concentration of free calcium ions in the blood plasma "ionized calcium". Siggaard-Andersen O, ed. *Blood pH, carbon dioxide, oxygen and calcium-ion*. 1 ed. Copenhagen: Private Press, 1981: 163-90.
25. Toffaletti J, Ernst P, Hunt P, Abrams B. Dry electrolyte-balanced heparinized syringes evaluated for determining ionized calcium and other electrolytes in whole blood. *Clin Chem* 1991; 37,10: 1730-33.
26. Fincher R, Strong J, Jackson J. Accuracy of measurements of hemoglobin and potassium in blood samples from peripheral catheters. *Am J Crit Care* 1998; 7,6: 439-43.
27. Spandrio L. Pre-analysis variability: the effect of evaporation on the error of measurement, with special reference to sodium and potassium. *Quad Sclavo Diagn* 1984; 20, 1: 110-21.
28. Kost GJ, Nguyen TH, Tang Z. Whole-blood glucose and lactate. *Arch Pathol Lab* 2000; 124: 1128-34.
29. Young DS. Effects of preanalytical variables on clinical laboratory test. 2nd ed. Washington, DC: AACC Press, 1997.
30. Reference Manual for the ABL700 Series. Ed. I. Radiometer publication, Copenhagen. Radiometer Medical A/S, 2003.
31. Andersen O, Haugeaard SB, Jørgensen LT *et al.* Preanalytical handling of samples for measurement of plasma lactate in HIV patients. *Scand J Clin Lab Invest* 2003; 63, 449-54.
32. Astles R, Williams C, Sedor F. Stability of plasma lactate in vitro in the presence of antiglycolytic agents. *Clin Chem* 1994; 40, 7: 1327-30.
33. Calatayud O, Tenias JM. Effects of time, temperature and blood cell counts on levels of lactate in heparinized whole blood gas samples. *Scand J Clin Lab Invest* 2003; 63: 311-14.

Andre nyttige ressursar:

Video: Arterial Puncture for Blood Gas Analysis. Shelly P. Dev, M.D., Melinda D. Hillmer, M.D., BSc.Ph.m., and Mauricio Ferri, M.D. [February 3, 2011](https://www.nejm.org/doi/full/10.1056/NEJMvcm0803851#figure=vcm0803851_thumb111x111.jpg) *N Engl J Med* 2011; 364:e7 DOI: 10.1056/NEJMvcm0803851: https://www.nejm.org/doi/full/10.1056/NEJMvcm0803851#figure=vcm0803851_thumb111x111.jpg

<https://www.uptodate.com/contents/arterial-blood-gases>

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Innsetting av venekateter. Illustrasjon © K. Toverud